

Alginate-Methacrylate Photoresins: Source, M/G Composition, Methacrylation Route, Covalent-Ionic Dual Networks, Hydration, Degradation and Qualification

Juan Seguro

3Dresyns by Resyner Technologies S.L.

Registered tax address: General Etxague 15B, 20003 Donostia / San Sebastián, Gipuzkoa, Spain

VAT number: ESB24827511

Correspondence: info@3Dresyns.com · Website: www.3Dresyns.com

ER-125 | Engineering Review | Version 0.3 | 31 August 2026

Statement of Significance

Alginate methacrylate is not one material. Algal source, molecular-weight distribution and mannuronate/guluronate composition define the incoming alginate; the attached methacrylate group and measured substitution define a photoactive derivative; calcium can create a second, reversible ionic network; and washing, ion exchange and hydration create later article states. ER-125 defines a qualification chain that prevents ALMA, methacrylated alginate and calcium-treated constructs from being pooled by name.

Keywords: alginate; alginate methacrylate; ALMA; mannuronate; guluronate; degree of substitution; calcium; ionic crosslinking; covalent network; dual network; hydration; degradation; qualification.

Abstract

Primary studies show that photocrosslinkable alginates are route-specific macromers, not interchangeable grades. Carbodiimide coupling of 2-aminoethyl methacrylate enabled substitution-dependent swelling, modulus and degradation; glycidyl methacrylate, 2-aminoethyl methacrylate and methacrylic anhydride produced derivatives with similar measured substitution but different compressive failure behaviour; functionalisation time and initiator concentration changed a separate photocurable bioink system; oxidation before methacrylation created another degradation coordinate; and light-patterned ionic alginate demonstrated that divalent-cation-bound networks can be fabricated without covalent methacrylate chemistry. Qualification must therefore freeze alginate source, molecular-weight distribution, M/G composition, modification route, measured substitution, oxidation, purification, concentration, ions, initiator, absorber, dose, working curve, conversion, wash, hydration, dimensions, mechanics, degradation protocol and model-specific biology. Calcium treatment is not proof of a covalent dual network, exposure time is not conversion, and one degradation profile is not a universal resorption rate.

1. Scope and boundaries

ER-125 owns material-specific identity and qualification of alginate methacrylate photoresins and explicitly declared covalent-ionic alginate networks. It includes alginate functionalised through glycidyl methacrylate, 2-aminoethyl methacrylate or methacrylic-anhydride routes, and oxidised-methacrylated alginate when oxidation is separately controlled. It does not own unmodified alginate generally, every calcium alginate gel, generic hydrogel architecture, universal biodegradation, tissue regeneration or clinical suitability.

ER-06 retains biological response; ER-10 geometric metrology; ER-15 degradation and mechanical evolution; ER-27 conversion; ER-33 washing and drying; ER-67 generic double-network architecture; ER-86 extractables; ER-100 corpus architecture; and ER-121 to ER-124 the neighbouring silk, collagen, chitosan and hyaluronic-acid derivative families [C1,C2,C3,C4,C5,C6,C7,C8,C9,C10,C11,C12].

M/G composition is not molecular weight. Modification feed is not measured substitution. Glycidyl-methacrylate, aminoethyl-methacrylate and methacrylic-anhydride derivatives are not chemically identical. Calcium-crosslinked native alginate is not ALMA, and an ALMA construct exposed to calcium is not a demonstrated dual network unless both covalent and ionic contributions are verified.

2. Controlled alginate-photoresin state chain

State	Minimum identity	Failure if omitted
Incoming alginate	Species/source route, salt form, molecular-weight distribution, M/G and block composition, purity, endotoxin/bioburden, storage	The word alginate replaces polymer identity
Photoactive derivative	Reagent and bond route, feed, pH, temperature, time, oxidation, purification and measured substitution	Feed or ALMA acronym is treated as chemistry
Vat or bioink	Derivative concentration, water/buffer, ions, initiator, absorber, cells/additives, rheology, temperature and age	Formulation and ion state disappear
Released article	Dose, conversion, wash, calcium history, hydration medium/time, dimensions, mechanics, residuals and degradation state	Green, ionic and hydrated articles are pooled

Table 1. Minimum state chain for alginate-derived photocurable networks.

This material/article chain is separate from process and production-site qualification under ISO/ASTM 52920:2023 [S2].

3. Primary evidence and controlled interpretation

3.1 Substitution-dependent covalent networks

Jeon et al. coupled 2-aminoethyl methacrylate to sodium alginate through carbodiimide/NHS chemistry and photocrosslinked the resulting macromers with ultraviolet light and 0.05% photoinitiator [1]. Swelling, elastic modulus and degradation changed with degree of methacrylation. The result establishes substitution as a controlled coordinate within that derivative and protocol; it does not make initiator percentage, ultraviolet exposure or degradation transferable to another alginate source or attachment chemistry.

3.2 Similar substitution does not erase attachment chemistry

Araiza-Verduzco et al. prepared alginate-glycidyl methacrylate, alginate-2-aminoethyl methacrylate and alginate-methacrylic anhydride derivatives from a 120-180 kDa alginate with a reported M/G ratio of 1.42 [2]. Their NMR calculation gave 31.39%, 28.03% and 32.35% substitution, respectively. At 2.5% w/v macromer, 10% w/w I2959 relative to macromer and 365 nm irradiation for 1 h, Young's moduli were about 15 kPa without a significant difference, while compressive strength and toughness followed a route-dependent order. Similar substitution and elastic response therefore did not prove chemical or failure equivalence.

3.3 Functionalisation time and initiator are separate variables

Mishbak et al. varied alginate functionalisation time and prepared photocurable systems with 0.5-1.5% w/v photoinitiator, followed by ultraviolet exposure at 8 mW/cm² for 8 min [3]. They evaluated rheology, mechanics, morphology, swelling and degradation. Those settings define the tested system; functionalisation time is not a substitute for measured substitution, and an eight-minute exposure is not a conversion result or a DLP working curve.

3.4 Ionic light-patterning is a distinct chemistry

Valentin et al. used stereolithographic patterning with photoacid generators and divalent cation salts to form ionically crosslinked sodium-alginate structures with reversible patterning and controlled degradation [4]. That study is important as a boundary: spatially fabricated divalent-cation-crosslinked alginate is not evidence for methacrylate conversion. Conversely, adding calcium after ALMA exposure does not establish two load-bearing networks without matched covalent-only, ionic-only and combined controls.

3.5 Oxidation creates another independent coordinate

Jeon et al. oxidised alginate before methacrylation and showed that oxidation level, polymer concentration and time altered swelling, modulus evolution and degradation of the resulting photocrosslinkable hydrogels [5]. Oxidised-methacrylated alginate must therefore be identified separately from unoxidised ALMA. Oxidation degree, actual methacrylation, molecular-weight change and chain-scission history cannot be collapsed into one nominal degradation grade.

4. Photochemistry, ions and process state

Alginate concentration, molecular weight, M/G/block distribution, buffer, salt and temperature alter viscosity and chain conformation. Substitution route, substitution level, initiator, oxygen, cells and absorber alter radical and optical response. Working curves must be measured on the exact vat state; exposure time alone is not conversion, and cure depth does not establish x/y/z fidelity [S1,S3,S4,C2,C4].

For a claimed covalent-ionic dual network, the sequence is part of the material identity: ionic pre-gel then light, light then calcium, or simultaneous exposure can produce different transport, shrinkage, surface and stress states. Calcium species, activity/concentration, counter-ion, bath volume, temperature, time and exchange history must be reported. Chelator-triggered or ion-exchange softening is not covalent degradation.

5. Hydration, degradation and biological boundaries

Evidence block	What it can show	What it cannot establish alone
Hydration/equilibrium	Medium- and time-specific mass, dimensions, swelling and mechanics	Conversion, in-use dimensions or degradation lifetime
Calcium/ion exposure	Ion-, bath-, time- and geometry-specific ionic contribution	Covalent conversion or a universal dual network
Hydrolysis/oxidation	Route-, pH-, medium-, temperature- and time-specific chain/network evolution	Universal in-vivo resorption
Cell/animal model	Response of the exact released construct in the declared model	Universal biocompatibility, efficacy or approval

Table 2. Evidence blocks that must remain separate for alginate photo-networks.

ISO 10993-12:2021/Amd 1:2025 and ISO 10993-18:2020/Amd 1:2022 anchor sample preparation and chemical characterisation [S5,S6]. ISO 10993-1:2025 frames biological-safety evaluation within risk management [S7]. Where an alginate hybrid contains an animal-derived component, ISO 22442-1:2020 frames source-risk management [S8]. Research-cell or animal observations do not constitute approval.

6. Release matrix

Gate	Minimum evidence	Claim withheld
Alginate identity	Source/species, salt form, molecular-weight distribution, M/G-block composition, purity and storage	Generic alginate equivalence
Functionalisation	Bond route, reaction conditions, measured substitution method/result, oxidation and residual-reagent control	Feed ratio or ALMA acronym as identity
Network attribution	Covalent-only, ionic-only and combined controls with sequence and ion history	Unverified dual-network reinforcement
Photochemistry	Spectrum, irradiance, working curve, conversion, gel fraction and spatial state	Exposure time as complete cure
Article state	Wash, ion exchange, hydration, x/y/z dimensions, mechanics and residuals	Green values as wet-duty properties
Biological duty	Defined degradation, cell and animal protocols on the released state	Universal safety, resorption or efficacy

Table 3. Minimum release architecture for alginate-derived photocurable constructs.

7. Failure-mode interpretation

- A nominally identical ALMA becomes much more viscous: verify source, molecular-weight distribution, M/G-block composition, concentration, salt, temperature and ageing before changing exposure.
- A higher methacrylation feed does not increase modulus: measure incorporated substitution and conversion; do not treat reagent feed as network density.
- Calcium treatment increases early modulus but the effect decays in buffer: map ion exchange and retain ionic attribution rather than relabelling it covalent degradation.
- The green construct is accurate but swells anisotropically: report x/y/z after a defined wash, ion treatment and equilibration.
- Oxidised ALMA loses mass faster: separate oxidation, molecular-weight change, substitution, conversion and chain-scission history before assigning one mechanism.
- Cells remain viable after printing: retain the exact derivative, light, ion, wash and culture state; do not convert one model into universal biocompatibility.

8. Minimum workflow and research priorities

- Freeze alginate species/source, salt form, molecular-weight distribution, M/G and block composition, purity, endotoxin/bioburden and storage.
- Report attachment reagent and bond route, feed, pH, temperature, time, oxidation, purification and measured substitution with its calculation basis.
- Map viscosity, optical response and working curves across temperature, salt state and the complete build-duration window.
- Measure conversion, gel fraction, residuals and retained alginate on matched green, washed, ion-treated and hydrated states.
- Use covalent-only, ionic-only and combined controls to verify a dual network and its sequence dependence.
- Track mass and x/y/z dimensions through washing, calcium exposure, equilibration, ion exchange and degradation.
- Report degradation medium, pH, temperature, replenishment, geometry and time rather than one generic rate.
- Separate model-specific viability or matrix observations from safety, efficacy and regulatory conclusions.

9. Conclusions

Alginate-derived photocurable hydrogels are a family of controlled polymer, derivative, network and article states. Source, molecular weight and M/G-block composition define the incoming alginate; attachment route and measured substitution define photoactivity; light defines a covalent network state; calcium and ion exchange can define a second, reversible state; and washing, hydration and degradation create later articles.

Primary evidence supports substitution-dependent tuning, route-dependent failure behaviour, functionalisation-sensitive photoinks, oxidised derivatives and light-patterned ionic alginate. It also shows why ALMA, exposure time, calcium treatment and degradation cannot be used as complete identities.

Release therefore requires traceability from alginate source through composition, molecular weight, functionalisation, oxidation, purification, formulation, ions, photochemistry, washing, hydration, dimensions, residuals, mechanics, degradation and model-specific biology. No single print, modulus, swelling ratio or cell result closes that chain.

Limitations

This paper is a qualitative engineering review, not a systematic meta-analysis. The cited studies use different alginate sources, molecular weights, M/G compositions, modification routes, substitution methods, oxidation states, concentrations, initiators, wavelengths, exposures, ions, geometries, hydration states, degradation protocols and biological models. Numerical results are not pooled or ranked. No proprietary formulation, universal DLP window, medical-device qualification or clinical claim is provided.

Author contributions

Juan Seguro: conceptualisation, evidence synthesis, boundary definition, engineering interpretation, manuscript preparation and final author verification.

Competing interests

The author is founder and CEO of 3Dresyns by Resyner Technologies S.L., which develops photopolymer materials. This relationship is declared for transparency. No commercial product claim is evaluated in this paper.

Data availability

No new experimental dataset was generated. Scientific statements derive from the cited primary literature and official standards. The review package includes QA and source-boundary records for independent review.

Company information

3Dresyns by Resyner Technologies S.L. · General Etxague 15B, 20003 Donostia / San Sebastián, Gipuzkoa, Spain · VAT ESB24827511 · info@3Dresyns.com · www.3Dresyns.com

References - independent scientific evidence

- [1] Jeon, O., Bouhadir, K.H., Mansour, J.M. and Alsberg, E. Photocrosslinked Alginate Hydrogels with Tunable Biodegradation Rates and Mechanical Properties. *Biomaterials* 30(14) (2009) 2724-2734. DOI: 10.1016/j.biomaterials.2009.01.034. PMID: 19201462. <https://doi.org/10.1016/j.biomaterials.2009.01.034>
- [2] Araiza-Verduzco, F., Rodríguez-Velázquez, E., Cruz, H., Rivero, I.A., Acosta-Martínez, D.R., Pina-Luis, G. and Alatorre-Meda, M. Photocrosslinked Alginate-Methacrylate Hydrogels with Modulable Mechanical Properties: Effect of the Molecular Conformation and Electron Density of the Methacrylate Reactive Group. *Materials* 13(3) (2020) 534. DOI: 10.3390/ma13030534. PMID: 31979162. <https://doi.org/10.3390/ma13030534>
- [3] Mishbak, H.H., Cooper, G. and Bartolo, P.J. Development and Characterization of a Photocurable Alginate Bioink for Three-Dimensional Bioprinting. *International Journal of Bioprinting* 5(2) (2019) 189. DOI: 10.18063/ijb.v5i2.189. PMID: 32954039. <https://doi.org/10.18063/ijb.v5i2.189>
- [4] Valentin, T.M., Leggett, S.E., Chen, P.-Y., Sodhi, J.K., Stephens, L.H., McClintock, H.D., Sim, J.Y. and Wong, I.Y. Stereolithographic Printing of Ionically-Crosslinked Alginate Hydrogels for Degradable Biomaterials and Microfluidics. *Lab on a Chip* 17(20) (2017) 3474-3488. DOI: 10.1039/C7LC00694B. PMID: 28906525. <https://doi.org/10.1039/C7LC00694B>
- [5] Jeon, O., Alt, D.S., Ahmed, S.M. and Alsberg, E. The Effect of Oxidation on the Degradation of Photocrosslinkable Alginate Hydrogels. *Biomaterials* 33(13) (2012) 3503-3514. DOI: 10.1016/j.biomaterials.2012.01.041. PMID: 22336294. <https://doi.org/10.1016/j.biomaterials.2012.01.041>

Standards and measurement references

- [S1] ISO/ASTM 52927:2024. Additive manufacturing - General principles - Main characteristics and corresponding test methods. Edition 1. <https://www.iso.org/standard/81802.html>

- [S2] ISO/ASTM 52920:2023. Additive manufacturing - Qualification principles - Requirements for industrial additive manufacturing processes and production sites. Edition 1. <https://www.iso.org/standard/76911.html>
- [S3] ISO/ASTM 52902:2023. Additive manufacturing - Test artefacts - Geometric capability assessment of additive manufacturing systems. Edition 2. <https://www.iso.org/standard/79683.html>
- [S4] ISO 1101:2017. Geometrical product specifications (GPS) - Geometrical tolerancing - Tolerances of form, orientation, location and run-out. Edition 4. <https://www.iso.org/standard/66777.html>
- [S5] ISO 10993-12:2021 with Amendment 1:2025. Biological evaluation of medical devices - Part 12: Sample preparation and reference materials. Edition 5. <https://www.iso.org/standard/75769.html>
- [S6] ISO 10993-18:2020 with Amendment 1:2022. Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process. Edition 2. <https://www.iso.org/standard/64750.html>
- [S7] ISO 10993-1:2025. Biological evaluation of medical devices - Part 1: Requirements and general principles for the evaluation of biological safety within a risk management process. Edition 6. <https://www.iso.org/standard/10993-1>
- [S8] ISO 22442-1:2020. Medical devices utilizing animal tissues and their derivatives - Part 1: Application of risk management. Edition 3; moved to stage 90.92 (to be revised); replacement ISO/AWI 22442-1 under development. <https://www.iso.org/standard/74280.html>

Related Engineering Review boundaries

- Qualification Constraints in Dental and Biomedical Applications. ER-06. Zenodo (2026). chemistry, Workflow Variables and
- [C2] Segurola, J. Metrology of Vat Photopolymerisation: Resolution, Trueness, Precision and Geometric Fidelity Across Scales. ER-10. Zenodo (2026).
- Qualification. ER-15. Zenodo (2026). al Evolution in Vat-Printed Tissue-Engineering Scaffolds: Coupled Design and
- Zenodo (2026). eterogeneity in Vat Photopolymerisation: From Exposure to Network State. ER-27.
- Zenodo (2026). g-Induced Extraction, Property Change and Dimensional Response in Vat-Printed Photopolymers. ER-33.
- Dissipation. ER-67. Zenodo (2026). gels in Vat Photopolymerisation: Microgel Reinforcement, Interpenetration and Energy
- and Cross-Study Comparability. ER-86. Zenodo (2026). opolymers: Specimen State, Extraction Design, Analytical Coverage
- ER-100. Zenodo (2026). ence Architecture Across Engineering Reviews ER-01 to ER-99.
- State and Qualification. ER-121. Zenodo (2026). , GMA Functionalisation, Aqueous DLP, Hydration, Mechanics, Biological
- Enzymatic Degradation and Qualification. ER-122. Zenodo (2026). Triple-Helix State, Low-Temperature DLP, Hydration,
- Hybrid PEG Networks, SLA/4D Printing and Qualification. ER-123. Zenodo (2026). n, Methacrylation, Acid-Neutral State,
- Hydration, Hyaluronidase Degradation and Qualification. ER-124. Zenodo (2026). Substitution, DLP Hybrid Networks,

Version note

v0.3 closes the standards-title and editorial findings and incorporates the independent-review corrections. The consolidated independent review dated 31 August 2026 returned APPROVE with the remaining primary-source verification closed from the enclosed full texts.

SHA-256