

Photocurable Pectin-Derived Hydrogels: Botanical Identity, Esterification, Methacrylation and Norbornene Routes, Ionic-Photo Dual Networks, Hydration, Degradation and Qualification

Juan Seguro

3Dresyns by Resyner Technologies S.L.

Registered tax address: General Etxague 15B, 20003 Donostia / San Sebastian, Gipuzkoa, Spain

VAT number: ESB24827511

Correspondence: info@3Dresyns.com | Website: www.3

ER-129 | Engineering Review | Version 0.3 | 31 August 202

Statement of Significance

Pectin is not a single reproducible macromer. Botanical source, extraction, galacturonic-acid content, molar-mass distribution, degree of methyl esterification, acetylation, salt and counterion state define the parent polymer; methacrylation or norbornene substitution defines a second identity; ionic and photo-crosslinking then create distinct network and article states.

Keywords: pectin; pectin methacrylate; norbornene-functionalised pectin; thiol-ene; dual network; bioink; hydration; qualification

Abstract

Photocurable pectin-derived systems include methacrylate-modified pectin combined with ionic pre-gelation, norbornene-functionalised pectin crosslinked with thiolated peptides, and hybrid dressings or bioinks in which pectin shares network ownership. Pereira et al. established a low-solids dual ionic-photo bioink and a separate thiol-norbornene platform; Giron-Hernandez et al. linked cocoa-pod extraction and esterification to pectin methacrylate; Bozer et al. studied a pectin/hyaluronic-acid hybrid. This review separates parent identity, derivatisation, formulation, network attribution, optical process, hydration, degradation and biological evidence.

1. Scope and boundaries

ER-129 owns material-specific identity and qualification of chemically photocurable pectin derivatives. It covers botanical source and extraction, pectin composition, methyl esterification, methacrylate or norbornene substitution, ionic pre-crosslinking, photochemistry, printing, washing, hydration, degradation, mechanics and biological evidence.

ER-06 retains biological-response methodology; ER-10 geometric metrology; ER-15 degradation and mechanical evolution; ER-27 conversion; ER-33 washing and drying; ER-67 generic dual-network architecture; ER-86 extractables; ER-100 corpus architecture; and ER-121 to ER-128 neighbouring natural-polymer derivative families. [C1,C2,C3,C4,C5,C6,C7,C8,C9,C10,C11,C12,C13,C14,C15,C16].

General AM characteristics and geometric test artefacts should follow the declared ISO/ASTM standards [S1,S3].

High- and low-methoxyl pectin are not interchangeable. Degree of methyl esterification is not degree of photoactive substitution. Ionic pre-gelation is not covalent conversion, and a pectin/nanocellulose extrusion ink with no pectin-bound photoactive group is not a photocurable pectin derivative. Methacrylate chain-growth and thiol-norbornene step-growth networks require separate controls.

2. Controlled pectin-photoresin state chain

State	Minimum identity	Failure if omitted
Incoming pectin	Botanical source/part, extraction and purification, lot, galacturonic-acid content, molar-mass distribution, methyl esterification, acetylation, ash, salt and storage	Pectin name replaces feedstock identity
Photoactive derivative	Methacrylate or norbornene route, feed basis, reaction, purification, measured substitution, residuals and storage	Feed ratio or derivative acronym is treated as identity
Ink or precursor	Pectin solids, counterions, co-polymers, peptides/cells, pH, ions, initiator, optical state, temperature and age	Formulation and network ownership disappear

State	Minimum identity	Failure if omitted
Article	Ionic/photo sequence, exposure, post-cure, wash, green/hydrated dimensions, conversion, soluble fraction, residuals, mechanics and degradation	Pre-gelled, covalent, washed and hydrated states are conflated

Table 1. Minimum state chain for photocurable pectin-derived networks.

This material and article chain is separate from process and production-site qualification under ISO/ASTM 52920:2023 [S2].

3. Primary evidence and controlled interpretation

3.1 Parent pectin identity begins before derivatisation

Giron-Hernandez et al. recovered pectin from cocoa pod husk by an optimised alkaline route at about 20% yield and reported about 47% degree of esterification before methacrylation [3]. Pereira et al. instead used purified low-methoxyl citrus pectin with 86% galacturonic-acid-unit content and 37% degree of methylation for norbornene functionalisation [2]. Botanical source, extraction yield, galacturonic-acid content and esterification therefore remain explicit lot coordinates, not generic pectin descriptors.

3.2 Methacrylated pectin enables a controlled ionic-photo sequence

Pereira et al. designed methacrylate-modified pectin bearing integrin-binding motifs and formed hydrogels by ionic gelation followed by UV photopolymerisation [1]. Ionic crosslinking adjusted the rheology of a 1.5 wt% solution so deposited filaments retained shape before post-print photo-crosslinking. The result is a pectin-based dual network: calcium or other ionic state, photoactive substitution, conversion and sequence must all be reported before printability or mechanics can be transferred.

3.3 Source-to-bioink traceability is experimentally achievable

Giron-Hernandez et al. confirmed conversion of extracted cocoa pectin to PECMA by FTIR-ATR and NMR, then combined CaCO₃ crosslinking with UV exposure [3]. UV changed porosity and improved compressive and rheological responses relative to the compared condition, while chondrocytes remained viable over the reported six-day culture. These outcomes belong to that waste source, extraction, esterification, derivatisation, ionic system, light history and cell model; they are not universal PECMA values.

3.4 Norbornene substitution creates a different reaction and degradation logic

Pereira et al. converted the low-methoxyl pectin to its tetrabutylammonium salt, reacted it with carbic anhydride in DMSO and measured 20.0 plus or minus 0.6% norbornene substitution by proton NMR [2]. A 2.5 wt% NorPEC formulation with 0.25 wt% VA-086 and 4-6 mM MMP-cleavable peptide was exposed at 365 nm and 7 mW/cm² for 20 s. Changing peptide crosslinker shifted swollen shear modulus from about 70 to 447 Pa and altered collagenase-triggered disintegration, demonstrating that linker identity and stoichiometry own degradability as well as network density.

3.5 Hybrid and non-photoactive pectin systems require explicit exclusion labels

Bozer et al. combined methacrylated pectin with methacrylated hyaluronic acid and oleuropein in a photocrosslinked wound-dressing study [4]; its wound-healing and cytocompatibility results describe a hybrid network, not pectin alone. Teixeira et al. bioprinted pectin, nanocellulose and lysozyme-nanofibril hydrogel bioinks [5]; unless a pectin-bound photoactive group and its conversion are reported, such a system is a pectin-containing composite and not evidence of a covalently photocurable pectin derivative. These classes stay outside neat-PECMA or NorPEC attribution unless matched component and network controls are supplied.

4. Photochemistry, ionic state and process control

Pectin source, galacturonic-acid fraction, molar-mass distribution, methyl esterification, acetylation, salt, counterion, pH, ionic strength and storage control solution and ionic-gel state. Photoactive substitution, initiator, wavelength, irradiance, absorber, exposure, oxygen management and temperature control covalent cure. Both coordinate sets must be measured on the production formulation over its use-life.

For dual-crosslinked PECMA, an ionic pre-gel can improve deposition yet obscure the contribution of photoconversion. For NorPEC, thiol:norbornene stoichiometry, pendant ligand content and degradable-linker concentration are network variables. Native pectin, ionic-only, photo-only, dual-network and hybrid-polymer controls are needed to assign rheology, conversion, stiffness, swelling and degradation.

5. Hydration, degradation and biological boundaries

Evidence block	What it can show	What it cannot establish alone
Parent identity	Method-specific source, composition, molar mass, methyl esterification and salt state	Photoactive substitution, conversion or article performance

Evidence block	What it can show	What it cannot establish alone
Ionic/process window	Formulation-, counterion-, printer- and geometry-specific deposition and cure behaviour	Transferable print settings or covalent conversion
Hydration/degradation	Medium-, ion-, enzyme- and time-specific mass, dimensions and mechanics	Universal degradation mechanism or service life
Biological model	Response of the exact construct or extract in the stated model	General safety, clinical efficacy or regulatory approval

Table 2. Evidence blocks that must remain separate for pectin-derived photoresins.

ISO 10993-12:2021/Amd 1:2025 and ISO 10993-18:2020/Amd 1:2022 anchor sample preparation and chemical characterisation [S4,S5]. ISO 10993-1:2025 frames biological safety within risk management [S6]. Cell viability in one construct, an extract assay or a scratch test does not establish universal biocompatibility, wound-healing efficacy or regulatory approval.

6. Release matrix

Gate	Minimum evidence	Claim withheld
Parent identity	Botanical source/part, extraction, lot, composition, molar mass, methyl esterification, acetylation, salt and storage	Generic pectin equivalence
Derivatisation	Reaction and purification record, measured photoactive substitution, residuals and stability	Feed or acronym as derivative identity
Process	Complete formulation, ionic sequence, printer/light hardware, wavelength, irradiance, exposure and geometry	Transferable deposition or cure settings
Network attribution	Conversion plus native, ionic-only, photo-only, dual and matched hybrid controls	Unverified ionic, covalent or co-polymer ownership
Article state	Green, post-cured, washed, hydrated and dried x/y/z, mass, residuals, mechanics and conditioning	Unqualified dimensions, mechanics or degradation
Biological boundary	Exact derivative, formulation, wash/extract, model, endpoint, duration and statistics	Universal biocompatibility or clinical suitability

Table 3. Minimum release architecture for photocurable pectin-derived constructs.

7. Failure-mode interpretation

- The nominally same pectin solution pre-gels faster: check botanical lot, molar-mass distribution, methyl esterification, free carboxylates, salt/counterion, calcium history, pH, ionic strength, temperature and age before changing exposure.
- More methacrylic anhydride gives no expected modulus gain: measure incorporated photoactive substitution and conversion; methyl esterification, reagent feed, grafted functionality, ionic junction density and covalent network density are distinct quantities.
- A deposited filament holds shape but loses form after washing: separate ionic pre-gel strength from covalent conversion, then quantify counterion exchange, soluble fraction, mass and x/y/z dimensions through the full wash and hydration sequence.
- A NorPEC gel is stiff but cells do not remodel it: verify thiol:norbornene balance, peptide identity, incorporation and accessible cleavage sequence; high crosslink density can slow enzyme-triggered network loss even when initial conversion is high.
- A hydrated article softens or swells across lots: condition to equilibrium in the intended medium and report ion composition, pH, time, temperature, mass, dimensions and mechanics; dry-state values do not qualify a hydrated pectin network.
- A cell-laden or wound-dressing model passes: retain the exact derivative, co-polymer, initiator, light dose, wash, cell source, exposure duration and endpoint; do not infer clinical safety, efficacy or regulatory approval.

8. Minimum workflow and research priorities

- Freeze botanical source, plant fraction, extraction and purification, lot, galacturonic-acid content, molar-mass distribution, methyl esterification, acetylation, ash, moisture, salt and storage before derivatisation.
- Report derivatisation reagent and feed basis, solvent, pH, time, temperature, atmosphere, purification, drying, measured methacrylate or norbornene substitution and calculation method.
- Map concentration, pH, ionic strength, calcium or other counterion, viscosity, yield stress, gel time, absorbance and working curves over temperature, ageing and build duration.
- Record printer and dispensing path, wavelength, irradiance, exposure, geometry, layer or filament conditions, delay between ionic and photo steps, post-cure and wash sequence.
- Measure photoactive conversion, insoluble and soluble fractions, residual initiator and derivatisation reagents, counterion content and x/y/z dimensions in green, washed, hydrated and dried states.

- Use native-pectin, ionic-only, photo-only, dual-network and matched hybrid controls to assign deposition, covalent cure, mechanical reinforcement and degradation.
- Test swelling, ion exchange, mass loss and mechanical retention in the intended medium; where degradable peptide linkers are claimed, include enzyme-free and relevant-enzyme controls.
- Separate extract cytotoxicity, direct contact, cell-laden printing, scratch assay, tissue model and animal observations from biological-safety, efficacy and regulatory conclusions.

9. Conclusions

Photocurable pectin materials are controlled botanical-feedstock, derivative, formulation, network, process and article states. Methyl esterification and photoactive substitution are independent, and ionic and covalent crosslinking require separate attribution.

Primary evidence supports dual ionic-photo PECMA bioinks, source-traceable cocoa PECMA, and peptide-crosslinked NorPEC with tunable hydrated mechanics and enzyme-responsive degradation. Hybrid PECMA/HA-MA and non-photoactive pectin composites are not neat-pectin evidence.

Release requires traceability through source, composition, esterification, derivatisation, ionic history, formulation, optical process, conversion, washing, soluble fraction, residuals, hydration, dimensions, mechanics, degradation and biological model.

Limitations

This paper is a qualitative engineering review, not a systematic meta-analysis. The cited studies use different botanical sources, extraction routes, esterification and substitution definitions, formulations, crosslinking mechanisms, optical processes, article states and biological models. Numerical values are study-specific evidence, not universal design values.

Author contributions

Juan Seguro: conceptualisation, evidence synthesis, boundary definition, engineering interpretation, manuscript preparation and final author verification.

Competing interests

The author is founder and CEO of 3Dresyns by Resyner Technologies S.L. No comparative product recommendation is made.

Data availability

No new experimental dataset was generated. Scientific statements derive from the cited primary literature and official standards. The review package includes QA and source-boundary controls.

Company information

3Dresyns by Resyner Technologies S.L. | General Etxague 15B, 20003 Donostia / San Sebastian, Gipuzkoa, Spain | VAT ESB24827511 | info@3Dresyns.com | www.3Dresyns.com

References - independent scientific evidence

- [1] Pereira, R. F., Sousa, A., Barrias, C. C., Bartolo, P. J. and Granja, P. L. A Single-Component Hydrogel Bioink for Bioprinting of Bioengineered 3D Constructs for Dermal Tissue Engineering. *Materials Horizons* 5(6) (2018) 1100-1111. DOI: 10.1039/C8MH00525G. <https://doi.org/10.1039/C8MH00525G>
- [2] Pereira, R. F., Barrias, C. C., Bartolo, P. J. and Granja, P. L. Cell-Instructive Pectin Hydrogels Crosslinked via Thiol-Norbornene Photo-Click Chemistry for Skin Tissue Engineering. *Acta Biomaterialia* 66 (2018) 282-293. DOI: 10.1016/j.actbio.2017.11.016. <https://doi.org/10.1016/j.actbio.2017.11.016>
- [3] Giron-Hernandez, J., Tombe, A., Chemban Koyilot, M., Salas-Calderon, K. T., Charlton, A., Wills, C. and Gentile, P. From Cocoa Waste to Sustainable Bioink: Valorising Pectin for Circular Economy-Driven Tissue Engineering. *European Polymer Journal* 210 (2024) 112967. DOI: 10.1016/j.eurpolymj.2024.112967. <https://doi.org/10.1016/j.eurpolymj.2024.112967>
- [4] Moran Bozer, B., Ozkahraman, B. and Mert, H. Photocrosslinked Methacrylated Pectin and Methacrylated Hyaluronic Acid Wound Dressing Loaded with Oleuropein as Bioactive Agent. *International Journal of Polymeric Materials and Polymeric Biomaterials* 73(11) (2024). DOI: 10.1080/00914037.2023.2235875. <https://doi.org/10.1080/00914037.2023.2235875>
- [5] Teixeira, M. C., Lameirinhas, N. S., Carvalho, J. P. F., Luis, J., Oliveira, H., Oliveira, J. M., Silvestre, A. J. D., Vilela, C. and Freire, C. S. R. Biobased Hydrogel Bioinks of Pectin, Nanocellulose and Lysozyme Nanofibrils for the Bioprinting of A375 Melanoma Cell-Laden 3D In Vitro Platforms. *International Journal of Biological Macromolecules* 282(Pt 5) (2024) 136958. DOI: 10.1016/j.ijbiomac.2024.136958. <https://doi.org/10.1016/j.ijbiomac.2024.136958>

Standards and measurement references

- [S1] ISO/ASTM 52927:2024. Additive manufacturing - General principles - Main characteristics and corresponding test methods. Edition 1. <https://www.iso.org/standard/81802.html>
- [S2] ISO/ASTM 52920:2023. Additive manufacturing - Qualification principles - Requirements for industrial additive manufacturing processes and production sites. Edition 1. <https://www.iso.org/standard/76911.html>
- [S3] ISO/ASTM 52902:2023. Additive manufacturing - Test artefacts - Geometric capability assessment of additive manufacturing systems. Edition 2. <https://www.iso.org/standard/79683.html>
- [S4] ISO 10993-12:2021 with Amendment 1:2025. Biological evaluation of medical devices - Part 12: Sample preparation and reference materials. Edition 5. <https://www.iso.org/standard/88205.html>
- [S5] ISO 10993-18:2020 with Amendment 1:2022. Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process. Edition 2. <https://www.iso.org/standard/64750.html>
- [S6] ISO 10993-1:2025. Biological evaluation of medical devices - Part 1: Requirements and general principles for the evaluation of biological safety within a risk management process. Edition 6. <https://www.iso.org/standard/10993-1>

Related Engineering Review boundaries

- [C1] Seguro, J. Biological Response in Photopolymer Additive Manufacturing: A Review of Formulation Chemistry, Workflow Variables and Qualification Constraints in Dental and Biomedical Applications. ER-06. Zenodo (2026).
- [C2] Seguro, J. Metrology of Vat Photopolymerisation: Resolution, Trueness, Precision and Geometric Fidelity Across Scales. ER-10. Zenodo (2026).
- [C3] Seguro, J. Transport, Degradation and Mechanical Evolution in Vat-Printed Tissue-Engineering Scaffolds: Coupled Design and Qualification. ER-15. Zenodo (2026).
- [C4] Seguro, J. Conversion Metrology and Spatial Cure Heterogeneity in Vat Photopolymerisation: From Exposure to Network State. ER-27. Zenodo (2026).
- [C5] Seguro, J. Washing- and Drying-Induced Extraction, Property Change and Dimensional Response in Vat-Printed Photopolymers. ER-33. Zenodo (2026).
- [C6] Seguro, J. Hierarchical and Double-Network Hydrogels in Vat Photopolymerisation: Microgel Reinforcement, Interpenetration and Energy Dissipation. ER-67. Zenodo (2026).
- [C7] Seguro, J. Analytical Extractables and Leachables from Vat-Printed Photopolymers: Specimen State, Extraction Design, Analytical Coverage and Cross-Study Comparability. ER-86. Zenodo (2026).
- [C8] Seguro, J. From Individual Reviews to a Qualification Corpus: Evidence Architecture Across Engineering Reviews ER-01 to ER-99. ER-100. Zenodo (2026).
- [C9] Seguro, J. Methacrylated Silk-Fibroin Photoresins: Protein Source, GMA Functionalisation, Aqueous DLP, Hydration, Mechanics, Biological State and Qualification. ER-121. Zenodo (2026).
- [C10] Seguro, J. Collagen Methacryloyl Photoresins: Collagen Source, Methacrylation, Triple-Helix State, Low-Temperature DLP, Hydration, Enzymatic Degradation and Qualification. ER-122. Zenodo (2026).
- [C11] Seguro, J. Chitosan-Derived Photocurable Hydrogels: Molecular Weight, Degree of Deacetylation, Methacrylation, Acid-Neutral State, Hybrid PEG Networks, SLA/4D Printing and Qualification. ER-123. Zenodo (2026).
- [C12] Seguro, J. Hyaluronic-Acid Methacrylate/Methacryloyl Photoresins: Source, Molecular Weight, Substitution, DLP Hybrid Networks, Hydration, Hyaluronidase Degradation and Qualification. ER-124. Zenodo (2026).
- [C13] Seguro, J. Alginate-Methacrylate Photoresins: Source, M/G Composition, Methacrylation Route, Covalent-Ionic Dual Networks, Hydration, Degradation and Qualification. ER-125. Zenodo (2026).
- [C14] Seguro, J. Dextran-Methacrylate Photoresins: Molecular Weight, Degree of Substitution, Solvent State, DLP Hybrid Networks, Washing, Degradation and Qualification. ER-126. Zenodo (2026).
- [C15] Seguro, J. Methacrylated-Starch Photoresins: Botanical Source, Amylose/Amylopectin Composition, Gelatinisation, Substitution, Aqueous DLP, Hybrid Bioinks and Qualification. ER-127. Zenodo (2026).
- [C16] Seguro, J. Photocurable Cellulose-Derived Photoresins: Cellulose and CMC Identity, Carboxymethyl and Methacrylate Substitution, Nanocrystal Functionalisation, Aqueous DLP/TPP, Hydration and Qualification. ER-128. Zenodo (2026).